



BIOSEN C_line

Menu structure of BIOSEN C_line




Info Sensor		
	Sensor - 1	Sensor - 2
Typ	Glucose	Lactate
Active	yes	yes
Calibration	502	374
Age	5	2

Info Counter	
	Running No.
Tray	1
C1	1
C2	1
Samples	1
STATs	1

Info Statistic		
	Today	Total
Tray	0	200
C1	0	18
C2	0	8
Samples	0	526
STATs	0	49

Info Device	
BIOSEN C_Line Clinic	
Model	5214-0121-0020
Ser. - No.	5214-04-0001
HMIC	1.02.0 / 400-002-018
PIC	1.02.0 / 400-002-015
	002-001



EKF- diagnostic GmbH
Ebendorfer Chaussee 3
Technologiepark Ostfalen
39179 Barleben/Magdeburg

Tel.: +49-(0)39203-785-0
Fax: +49-(0)39203-785-16
info@EKF-diagnostic.de
www.EKF-diagnostic.de

BIOSEN C-LINE

Demand-oriented
Glucose/Lactate
measurement